

DAVIDS EXPERIMENTAL LAB NOTEBOOK

CREATED: 09/07/2022


TABLE OF CONTENTS:

	<u>PAGES</u>
1. <u>Distillation of an unknown binary mixture</u>	2-10
a. Pre-Lab p.3-8	
b. Report/Data p.5-10	
2. <u>Recrystallization Methods and Melting Points</u>	11-24
a. Pre-Lab p.12-22	
b. Report/Data p.17-24	
3. <u>Extraction Methods</u>	25-34
a. Pre-Lab p.26-33	
b. Report/Data p.30-34	
4. <u>Optical Resolution of (+-) Mandelic Acid</u>	35-43
a. Pre-Lab p.36-42	
b. Report/Data p.43	
5. <u>Dehydration of Tetrahydrolinalool (Alkane->Alkene)</u>	44-52
a. Pre-Lab p.45-49	
b. Report/Data p.52	
6. <u>Williamson Ether Synthesis</u>	53-60
a. Pre-Lab p.53-60	
b. Report/Data p.60	
7. <u>End of Organic Chemistry 2022 Chapter</u>	61



[Return To Table of Contents](#)

DISTILLATION OF AN UNKNOWN BINARY MIXTURE



Experiment #1

Pre Lab: Distillation 2-ways

9/7/2022 4

- Title: Distillation Methods
- Objective: To use 2 techniques of distillation to purify, then identify an unknown liquid. The identification will be carried out using multiple diagnostic methods that will be cross-referenced with the properties of the possible chemicals.
- References: 1. ChemSurvival, Youtube. A Brief Introduction to Fractional Distillation. (2012, Jul).
2. Experiment 1, Distillation. Labflow.com. (N.D.)
3. Video - Simple Distillation. Labflow.com. (N.D.)
4. Video - Fractional Distillation. Labflow.com. (N.D.)

• Theory:

- Through distillation, one can effectively separate the azeotrope into its respective components. This is done by taking advantage of each miscible substances point of vaporization, and recondensing the pure vapor into a separate vessel... The purity can be enhanced using fractional distillation as the compound with the higher boiling point can get mixed in with the vapor of the desired compound. It is enhanced because fractional distillation essentially performs a simple distillation many times over through condensation and re-vaporization up the fractionating column.

- The boiling point of the distillate will be recorded over time and a graph will be created... by using refractometry along with the boiling point, we can deduce the identity of the compound by cross referencing it with the given table of possible compounds.

Pre-Lab: Distillation Methods

5

Theory cont.:

Chemical equations for rxns...

1. $(\text{Unknown} + \text{Mixture})_{(l)} \xrightarrow{\text{vaporization}} \text{unknown}_{(g)} + \text{Mixture}_{(g)}$
2. $\text{unknown}_{(g)} + \text{Mixture}_{(g)} \xrightarrow[\text{separation (reflux)}]{\text{distillation}} \text{unknown}_{(g)} + \text{Mixture}_{(g)}$
3. $\text{unknown}_{(g)} \xrightarrow{\text{condensation}} \text{unknown}_{(l)}$

Simple Distillation





Figure #1 - simple distillation

Fractional Distillation

This method uses the same set up as above, but adds a Fractionating column to the joint marked by a \star . This piece of equipment condenses the higher BP substance continuously and therefore further purifies the target distillate...

glass beads
increases to
surface area

Figure #2
Fractionating
Column.



Pre Lab: Distillation Methods

theory cont:

Calculations for Theoretical Yield:

- The theoretical yield can not be calculated beforehand due to the proportion of substance A & B being unknown... However, once they are known, ~~I can~~ ^{DEM} The following formula will be used

$$\text{Actual mL(A) given} \times \frac{g(A)}{\text{mL(A)}} \times \frac{\text{mol(A)}}{g(A)} = \text{Theo yield(A)}$$

~~% yield = $\frac{\text{Theo yield}}{\text{Measured Yield}}$~~

$$\% \text{Error} = \frac{\text{Theo yield} - \text{actual yield}}{\text{Theo yield}} \times 100\%$$

Tables of Reagents & Products

checkbox for elim.

Compound Name	Structure	Color	B. P. (°C)	n_D^{25}	Hazards info
<input checked="" type="checkbox"/> n-pentane		Colorless, clear	36	1.3547	mild irritant... Acute CNS depressant Acute sensory/motor nerve damage
<input checked="" type="checkbox"/> Acetone		Colorless, clear, odorless	56	1.3561	highly flammable serious eye irritation chronic exposure reorganizing
<input checked="" type="checkbox"/> Methanol		Colorless, clear	65	1.3288	highly flammable (15mg) causes organ damage toxic if swallowed
<input checked="" type="checkbox"/> n-Hexane		Colorless, clear	69	1.3723	highly flammable May be fatal if swallowed irritant
<input checked="" type="checkbox"/> Ethanol		clear, colorless	78	1.3610	Highly flammable causes serious eye irritation
<input checked="" type="checkbox"/> Cyclohexane		colorless, clear	81	1.4246	Highly flammable May be fatal if swallowed causes irritation
<input checked="" type="checkbox"/> 2-Methyl-2-Propanol		Colorless, clear	82	1.3878	highly flammable harmful by inhalation irritant
<input checked="" type="checkbox"/> Water		colorless, clear	100	1.3330	Non-Toxic You could drown if you inhale enough.

Table continued on next page

Pre-Lab: Distillation Methods

Checkboxes #1
for elimination

Table cont.

Compound Name	Structure	Color	B.P. (°C)	n_D^{25}	Hazards Info
<input checked="" type="checkbox"/> Methylcyclohexane	<chem>CC1CCCCC1</chem>	colorless clear	101	1.4231	Flammable, skin irrit., toxic to aquatic org. spp.
<input checked="" type="checkbox"/> 2-Methyl-1-propanol	<chem>CC(C)CO</chem>	colorless clear aromath. (A)	108	1.3937 1.3936	causes serious eye damage, flammable, skin irritant, may cause dizziness
<input checked="" type="checkbox"/> Toluene	<chem>Cc1ccccc1</chem>	colorless clear A	110	1.4941	highly flammable, may be fatal if swallowed, can damage fetus, may cause dizziness
<input checked="" type="checkbox"/> 1-Butanol	<chem>CCCCO</chem>	colorless clear Alc-like	117	1.3993	Flammable, skin irritant, may cause dizziness, harmful if swallowed
<input checked="" type="checkbox"/> m-Xylene	<chem>Cc1cc(C)ccc1</chem>	colorless clear	139	1.4941	highly flammable, serious eye irrit., may cause dizziness, chronic organ dmg.


Table #1 - Chemical Data table (possible chem)

References:

1. "Experiment 1" (2022), LabFlow.com, Storage-cbn, labflow.com/data/fileedit/0f/c0/0fc0c29d3f9#4409a860f928d
2. National Institute of Standards and Technology (NIST), N/D, webbook.nist.gov
3. Fisher Scientific, N/D, Fishersci.com.

Graph Ex. of Simple Dist. vs. Fractional Dist

50ml total
a 50/50 mix
of Acetone
& Cyclohexane
Acetone = 56°C
Cyc = 81°C



PART I Pre-Lab: Distillation Methods

Notes ✓ box Procedure Simple Distillation (we for Report) Data Simple Distillation

Binary Mixture!
 Obtain a 50 mL sample from Dr. G... It will be 2 of the chems in table in any proportion

Obtained 51.5 mL of unknown #6...
 - has an alcoholic odor slightly fruity...

Water in is the low side of condenser
 Ice bath?
 set up a simple distillation refer to p. 5 of ~~pre~~ lab. book ^{DEM}
 - use a 100 mL round bot. flask
 as a pot flask & a 10 mL graduated cylinder as a receiver.
 - use wood blocks to elevate receiver to a proper height.
 - be sure to grease joints!

- Set up apparatus...
 approved by Dr. G
 VARIAC was powered on 50 to begin boiling azeotrope... 120V was set.
 T shot up to 60°C
 - recorded 28 drops/min

Simple Dist

T ₉	mL
60.0	0
61.5	1
62.5	2
63.5	3
64.5	4
64.5	5
65.0	6
65.0	7
66.0	8
66.5	9
67.0	10
67.0	11
67.5	12
68.5	13
69.5	14
71.0	15
72.5	16
74.5	17
76.5	18
78.5	19
80.0	20
82.5	21
85.0	22
88.5	23
92.5	24
95.0	25
98.5	26
100.0	27
101.5	28
103.0	29
103.5	30
104.0	31
105.0	32
105.0	33
105.0	34
105.0	35
105.0	36

Transfer unknown liquid (50 mL) to the pot flask...
 - use a long-stemmed funnel
 - don't forget boiling chips
 Have Dr. G give the O.K.

- Cranked up to 60 .5 mL/min
 - cranked up to 65
 - 60 drops/min
 - Distillation complete
 46.75 mL of azeotrope left.

Use a heating mantle w/ the VARIAC to slowly bring the solution to a boil.
 - set VARIAC to about 40 and adjust as needed.

Don't forget to record every 1 mL of distillate.

Boil until about 2 mL of solution is left... Do not boil dry!
 - be sure to keep distillate rate @ 1-2 mL/min

Don't forget to change vials DEM as B.P. rises and collect the fractional dist later...
LATER

I repeat do not let to progress...
 Study SDS's

Table 1
 Create a data table

Pre-Lab ^{rem} Distillation Methods

9

Notes
✓ done

Procedure Simple Distillation cont.

- After finishing, take off the mantle and unplugging everything...
- Transfer the distillate back to the 100mL boiling (RB) flask.
 - cork it and store safely
 - graph data!

Clean all glass-ware!

Data simple Distillation cont.

- Side note: I almost distilled to dryness! there was only about .5mL left...

- Stop after ~5 constant temps next time and be sure to check under Al foil...

T°C	ml
105.0	37
105.5	38
105.5	39
105.5	40
105.5	41
105.5	42
105.5	43
105.5	44
105.5	45
109.5	46
109.5	47
STOP	48
Dist	49
Dist	50
Dist	50

PART II Fractional Distillation

of a binary mixture

- Set up a fractional dist. apparatus (p. 5 of lab book).
 - exactly the same as simple w/ a Vigreux coln.
 - Obtain 2 dry 25mL E. Flasks for pure samples...
 - maybe a third for the mix
 - Wrap column in glass wool or Al foil to drop heat.
- Repeat the procedure for simp. dist. with collections of A, AB, & B in separate flasks.

T°C @ 0mL = 56

T°C	ml	T°C	ml
57.0	1	109.0	26
57.5	2	105.0	27
57.5	3	105.0	28
58.0	4	105.0	29
58.5	5	105.0	30
60.0	6	105.0	31
60.5	7	105.0	32
60.5	8	105.0	33
61.0	9	105.0	34
61.0	10	105.0	35
61.0	11	105.0	36
61.5	12	105.0	37
61.5	13	105.0	38
61.5	14	105.0	39
62.0	15	105.0	40
62.0	16	105.0	41
60.5	17	105.0	42
61.0	18	105.0	43
61.0	19	105.0	44
63.0	20	105.0	45
103.5	21	105.0	46
105.0	22	105.0	47
105.5	23	105.0	48
105.5	24	105.0	49
104.5	25	105.0	50

F-Dist table for data

Tab/Graph Checklist

- Tab. S-Dist F-Dist RI-Index
- Graph S-Dist F-Dist
- OC vs. mL

Label and turn in samples w/ the reports.

- Determine Refractive index for A & B... using B.P. & RI, determine identity...
- Graph T vs. mL ^{rem}

Report: Distillation Methods

10

→ Summary Table

Distillates →	A (Acetone)	B (2-Methyl-1-propanol)
Observed B.P. (°C)	~ 61.0°C	105.0°C
R.I. (corrected) _{n_D^{25}}	1.35855	1.39310
Molar Mass (g/mol)	58.08 g/mol	74.12 g/mol

→ Critique of Experiment (ERROR)

Simple Distillation (Part I)

- As the data in the graph shows (and the table), the separation of A and B is not very clean and leaves a lot of room for A to contaminate B and vice versa... This is due to the vapors of the two compounds being able to "hitch a ride on each other."

- During experimentation, ~~50~~ 51.5^{DEM} mL of azeotrope was obtained and due to the small amount of vapor loss and pouring, there was only about 46.5 mL of solution at conclusion...

- The azeotrope was almost distilled to dryness and got a light brown coloration... This could have an effect on R.I. due to it being added to the main batch...

Fractional Dist. Error →

Report: Distillation Methods cont.

11

Fractional Distillation (Part II) ERROR

- loss of product due to residue in glassware and pouring
- Product was also lost (sacrificed) due to a pause in collection during the phase change period. (6.0 mL)
- Too much of the azeotrope was disposed of at the end (7.5 mL)... At a glance, the O.K. to remove from heat was given by instructor, but it was the fault of the experimenter to not use their judgment after observing the large quantity of solution remaining...
- The water for the condenser turned itself off and the lack of flow was not noticed for roughly 5-10 min.

X. David Miller

DATE ★ 09/20/2022.

★ Data explained in the excel data sheet.

Conclusions

By implementing knowledge of phase changes of chemicals into the process of distillation, Compound A & B were able to be successfully separated using

fractional distillation. To further increase the accuracy of identification, the boiling points were compared to the textbook values as well as refractive index. By using these values, A is Acetone and B is 2-methyl-1-propanol as shown on the table of reagents and the excel data sheet. This conclusion can be proven to have even more accuracy when referencing the values of purity through analysis of the liquids.

[Return To Table of Contents](#)

RECRYSTALLIZATION METHODS AND MELTING POINTS



Pre-Lab: Recrystallization Methods

12

★ Objective:

PI: To screen reaction solvents by testing the solubility of a number of known compounds.

PII: To carry out the crystallization of a known compound using the screening method from PI as a prerequisite to the experiment. A yellow impurity will also be purged from solution using Norit to yield a pure product.

PIII: To carry out the recrystallization of an unknown compound by using a "good" solvent followed by a "bad" solvent to optimize purity. Melting point analysis and infrared spectroscopy will be used to meet the objective of identifying the unknown.

★ References:

1. "Experiment 2" (2022). LabFlow.com

Storage - cdn.labflow.com/data/files/17/4c/174c69ac0693022

2. Fisherscientific, N/D. Fishersci.com

3. d21.wcupa.edu/d21/home/3444418

- Solubility table

Theory

Pre-Lab: Recrystallization Methods

13

* Theory

• Techniques:

- PI - Recrystallization is a technique used to remove impurities from a crystal sample to obtain the desired pure product. This is done by taking advantage of the crystal's solubility by testing it with different solvents (polar \rightarrow NE). A proper solvent is chosen by performing a screening of the solvents. This is performed by figuring out what solvent the crystal is insoluble with at room temperature and completely soluble with at hot/boiling temperature. After the solvent is chosen, the solution is heated up and the impurities are removed using hot gravity filtration through a fluted filter paper (heat is used to prevent crystallization of desired product). Activated charcoal (Norit) can also be used as a purifier (with an expense to final yield)...
- PII - Recrystallization can also be done using a method of a "good" and "poor" solvent with regards to the crystal. The "good" solvent (S @ R.T. & Hot) dissolves the crystal and is heated to boil. The "bad" solvent (I @ R.T. & Hot) is then added to the filtrate until cloudiness forms... (this forces the crystal out of solution). For all procedures, after the crystals form, they are extracted using a Büchner funnel in a vacuum filtration apparatus. This is done to separate out the crystals for the final yield.
- PIII - After the unknown has been recrystallized, one can determine its purity by using melting point ^{DETA} analysis. This is done by using a Mel-Temp to slowly raise T (C) until a melting point has been reached. The narrower the melting range, the purer the product. Cont. \longrightarrow


Pre-Lab: Recrystallization Methods

14

Theory cont... Techniques cont.

- Another method of analysis is infrared spectroscopy. This is a technique of firing infrared rays at a substance using different energy levels and analyzing the emission spectra of the resonating atomic characteristics. This method gives the experimenter information that can be used alongside the melting point to determine which of the possible crystals the unknown could be.

- Chemical Eqns. For rxns./side rxns.



Equation 1: ΔH_{diss} Dissolving of a polar Crystal solute

- Diagrams for Apparatus




Figure 1: Hot Gravity Filtration

Pre-Lab Cont.: Recrystallization Methods

15

Theory cont... Diagrams for apparatus cont.




Figure 2: Vacuum Filtration

- Calculations for Theoretical yield
 - T-Yield can't be accurately calculated due to not having data, however, the formula for % yield is as follows.
- $$\left[\frac{(T\text{-yield}) - (\text{actual yield})}{(T\text{-yield})} \right] \times 100\% = \% \text{ yield}$$

Table of Reagents: Solvents

Name	Structure	B.P. (°C)	MW (g/mol)	Density (g/ml)	Hazard info
Diethyl ether	<chem>CCOC</chem>	35°C	74 g/mol	0.71	<ul style="list-style-type: none"> • May damage fertility • Harmful if inhaled • Highly flammable • May cause organ damage
n-pentane	<chem>CCCCC</chem>	36°C	72	0.62	<ul style="list-style-type: none"> • Acute CNS depressant • Mild irritant
Dichloromethane	<chem>ClCCl</chem>	41°C	84	1.34	<ul style="list-style-type: none"> • Carcinogenic • May cause dizziness • May cause organ damage
Acetone	<chem>CC(=O)C</chem>	56°	58	0.79	<ul style="list-style-type: none"> • Flammable • May cause organ damage
Chloroform	<chem>ClC(Cl)Cl</chem>	61°	118	1.49	<ul style="list-style-type: none"> • Carcinogenic • Damages unborn child • Toxic if inhaled
Methanol	<chem>CH3OH</chem>	65°	32	0.79	<ul style="list-style-type: none"> • Flammable • Toxic if swallowed
n-Hexane	<chem>CCCCCC</chem>	69°	86	0.65	<ul style="list-style-type: none"> • Flammable • Harmful if swallowed

Pre-Lab: Recrystallization Methods

Table of Reagents: Solvents cont.

+ = "Good"
 - = "poor"
 H = "neither"

Name	Structure	B.P. (°C)	MW (g/mol)	Density (g/mL)	Hazard Info
+ Ethyl Acetate	<chem>CC(=O)OCC</chem>	77°	88	0.90	<ul style="list-style-type: none"> Flammable organ tox. Harmful if swall.
+ Carbon Tetrachloride	<chem>ClC(Cl)(Cl)Cl</chem>	77°	153	1.58	<ul style="list-style-type: none"> carcinogenic tox. if inhaled tox. to organs
+ Ethanol	<chem>CCO</chem>	78°	46	0.79	<ul style="list-style-type: none"> Flammable irritates eye
+ Benzene	<chem>c1ccccc1</chem>	80°	78	0.88	<ul style="list-style-type: none"> Flammable tox. if swallowed
+ Water	<chem>O</chem>	100°	18	1.00	N/A
+ Toluene	<chem>Cc1ccccc1</chem>	111°	92	0.86	<ul style="list-style-type: none"> Flammable tox. U.S. child tox. if swall. organ dmg.
+ Heptanes (Ligroin)	<chem>CCCCCCC</chem>	98°	87	0.68	<ul style="list-style-type: none"> Flammable CNS damage Fatal if swall. Organ dmg.
+ END	END OF SOLVENT TABLE				

Table 1.: Table of reagents: Solvents

Table of Reagents: Unknowns chart of Crystals

Name	Color	Structure	MP (°C)
------	-------	-----------	---------


Pre-Lab: Recrystallization Methods

Table of Reagents: Unknowns chart of Crystals

Name	Structure	MP(°C)	Color	Hazards
<input checked="" type="checkbox"/> Benzoic Acid	<chem>c1ccccc1C(=O)O</chem>	122°	Colorless	• Eye irritant
<input checked="" type="checkbox"/> Anthracene	<chem>c1ccc2cc3ccccc3cc2c1</chem>	216°	Colorless to pale yellow	• Flammable powder • eye irrit.
<input checked="" type="checkbox"/> Stearic Acid	<chem>CCCCCCCCCCCCCCCC(=O)O</chem>	70°	Colorless	• May form combustible dust
<input checked="" type="checkbox"/> Phthalic Acid	<chem>c1ccc(cc1C(=O)O)C(=O)O</chem>	206-208°	Colorless	• Causes serious eye damage
<input checked="" type="checkbox"/> Acetanilide	<chem>CC(=O)Nc1ccccc1</chem>	113-114°	Colorless	• tox. if swallowed
<input checked="" type="checkbox"/> Benzil	<chem>c1ccc(cc1)C(=O)C(=O)c2ccccc2</chem>	157° 95°	off-white yellow	• skin irrit • resp irrit • eye irrit
<input checked="" type="checkbox"/> Benzoin (AI)	<chem>c1ccc(cc1)C(=O)C(O)c2ccccc2</chem>	137°	off-white	N/A
<input checked="" type="checkbox"/> Diphenyl-acetic Acid	<chem>c1ccc(cc1)C(O)C(=O)c2ccccc2</chem>	146°	colorless	• causes skin irrit • causes eye irrit
<input checked="" type="checkbox"/> Triphenyl-methanol	<chem>c1ccc(cc1)C(O)(c2ccccc2)c3ccccc3</chem>	164-165°	colorless	• slightly flammable • eye irritant

Table 2: Table of Reagents - unknowns chart of Crystals

Procedure on next pg



Pre-Lab: Recrystallization Methods

18

Procedure

- Part I:
- Obtain 5 pipettes and 1 bulb
 - Label paper to keep track of pipettes.
 - Transfer 1 head of a small spatula (50mg) of powder to the test tubes and test against all solvents...
 - Follow elimination method learned in lab vids... Fill in table.
 - Heat all applicable solutions ^{test tubes} over a hot plate in a beaker with H₂O.
 - Recrystallize all suitable solutions in an ice bathed beaker. w/ ethanol
 - Clean + + and dump mixture in waste container.
 - scrub out w/ soap & H₂O

Data

I=insoluble SS=slightly soluble S=soluble

	Anthracene	Diphenylacetic Acid	Benzoin	Triphenyl Methanol	Name
H ₂ O	I	I	I	I	H ₂ O
EtOH	S	SS	SS	SS	EtOH
MeOH	S	SS	SS	SS	MeOH
EtAc	S	SS	SS	SS	EtAc
PhMe	S	SS	SS	SS	PhMe
PhOH	S	SS	SS	SS	PhOH
PhCl	S	SS	SS	SS	PhCl
PhI	S	SS	SS	SS	PhI
PhBr	S	SS	SS	SS	PhBr
PhNO ₂	S	SS	SS	SS	PhNO ₂
PhSO ₃ H	S	SS	SS	SS	PhSO ₃ H
PhCOOH	S	SS	SS	SS	PhCOOH
PhNH ₂	S	SS	SS	SS	PhNH ₂
PhNMe ₂	S	SS	SS	SS	PhNMe ₂
PhNMe	S	SS	SS	SS	PhNMe
PhNHMe	S	SS	SS	SS	PhNHMe
PhNHMe ₂	S	SS	SS	SS	PhNHMe ₂
PhNMe ₃	S	SS	SS	SS	PhNMe ₃
PhNHMe ₃	S	SS	SS	SS	PhNHMe ₃
PhNHMe ₂	S	SS	SS	SS	PhNHMe ₂
PhNHMe	S	SS	SS	SS	PhNHMe
PhNH	S	SS	SS	SS	PhNH
PhN	S	SS	SS	SS	PhN
Ph	S	SS	SS	SS	Ph
Temp	20°C	100°C	95°C	165°C	
Appearance	White/crystalline	colorless powder-like	yellow colored Seems to be non-p Benzoin floats in H ₂ O	colorless powdery	

Data (Expanded)

	Benzoin	Acetanilide	Temp
H ₂ O	I	S	20°C
EtOH	S	S	100°C
MeOH	S	S	95°C
EtAc	S	S	165°C
PhMe	S	S	
PhOH	S	S	
PhCl	S	S	
PhI	S	S	
PhBr	S	S	
PhNO ₂	S	S	
PhSO ₃ H	S	S	
PhCOOH	S	S	
PhNH ₂	S	S	
PhNMe ₂	S	S	
PhNMe	S	S	
PhNHMe	S	S	
PhNHMe ₂	S	S	
PhNMe ₃	S	S	
PhNHMe ₃	S	S	
PhNHMe ₂	S	S	
PhNHMe	S	S	
PhNH	S	S	
PhN	S	S	
Ph	S	S	

Table 3: Solubility Table

Pre-Lab: Recrystallization Methods

19

Procedure

Data

Part II

Obtain impure phthalic acid
(0.2% colored impurity)

Color: slightly yellow Texture: grainy, mud-like

weigh vial w/ cap
↓ transfer phthalic acid
weigh v + cap again

Weights

V+C+PA (g): 19.68g

↓
V+C (g): 17.64g

(g) Phthalic Acid = 2.04g

Use powder funnel to transfer
P-acid to E-Flask (125mL)
• return vial to D.G.

weigh roughly 300mg (10-15%)
Norit using weighing paper.

add 5mL DI water & boiling
stick to E-Flask... Boil it
add H₂O 1 drop at a time
until sample dissolves.

record total solvent needed

start solvent = 25.7 mL + 100% = 51.6 mL

Tot. solv. ^(H₂O) needed = 25.7 mL

↓ 2x

about = 51.8 mL

add 100% excess H₂O
remove from heat and add
Norit, swirl to mix

Boil for 5 min

Fold a fluted filter paper
and perform a hot gravity
filtration.

Matters: stainless funnel

• E-Flask • boiling stick



• fluted paper

• Ring stand

• a bit of solvent

(~1/3 full)

(before boiling)
Hot gravity solvent start = ~30 mL

1:50-1:55 pm

- Bubbles on outside have
formed w/ heat and solute
has dissolved...

- Added Norit @ 1:58
300 mg

- solution re-simmered
2:00 pm


- Remove soln. from heat

Pre-Lab: Recrystallization Methods

20

Part II
cont.

Procedure

- ⊗ Observe Filtrate color
 - ⊗ More Nitro needed?
 - ⊗ Color OK? (clear)
- ⊗ When colorless, cover flask w/ inverted beaker  and cool to R.T. slowly.
- ⊗ Cool in ice bath after crystal formation for highest crop.
- ⊗ Collect crystals using Vacuum Filtration.
 - Do not Filter over 2/3 full wet paper a bit w/ solvent before adding powder for maximum stick.
 - use ice cold solvent (H₂O) to rinse the last of the crystals out of E-flask
- ⊗ For a higher yield, evaporate filtrate to 2/3 volume by boiling and repeat process.
- ⊗ Transfer crystals to watch glass

- ⊗ Test MP, if ranges are identical, combine crops.

Data

- Performed Hot vac. filtration.
- Color of filtrate 1 = ~~##~~ Clear
- ⊗ Placed filtrate in ice bath

- Have crystals occurred?
- Solvent volume post filtration = ~75 mL

- No crystals formed so I ^{DEM} put back on heat to attempt to reduce down to 30 mL
- ⊗ They formed!

- Crystals are white
- performed vacuum filtration.
 - The vacuum busted a hole in my paper and I lost yield.
 - Some paper shreds maybe in yield.
- repeated recrystallization process using filtrate 2. ^{DEM}
- Watch glass = 47.28 g
- Dry weight crop 1 = 47.28 - 16.24 =
- crop 2 = N/A

- M.P. crop 1 = 203.1°C
- crop 2 = N/A
- % recovery crop 1 = 51.09%
- crop 2 = N/A
- crop 1 + 2 = ~~##~~ 60.4%

Pre-Lab: Recrystallization Methods

21

Procedure

Data

Part III

- Obtain unknown and note color/texture
- weigh unknown
- repeat solubility test from P I and determine "good" and "bad" solvents
- weigh unknown

after tests!

- add powder to 125 mL E-Flask
- add small amount of "good" solvent to bring to boil and cont. to add hot solvent until solute dissolves.

- add extra solvent to prevent crystals during hot grav. Filtration

- +50% for organic
- +100% for H₂O
- +25% for mixed solv. crystallization

add Norit if it is colored

- Repeat Hot gravity Filtration from P II
- don't forget the beaker



For MIXED

- add "poor" solvent drop by drop to hot filtrate. Keep all solvents hot. Do this until cloudy
- After crystallization occurs, use vacuum filtration

Unknown #33

Color: Sand colored (tan) Texture: Sugar-sized crystals (cone)

Vial + cap + X = 18.47g

Vial + cap (g) = 17.54g

(g) X = 0.92g

H₂O MeOH > EtOH Ligroin

RT SS	RT S	RT S	RT SS	RT I
S _B	S _B	S _B	S _B	SS _B

Total solvent added:

- performed a solubility test
- Ligroin was insoluble, leaving sandy bits behind
- MeOH and EtOH are viable for A
- MeOH > EtOH - H₂O is miscible

Total solvent added: (MeOH)

20 mL + 50% = 30 mL

Total Norit added = .28g

- heated solvent + solute Boil @ 2:30
- performed hot grav. Filtration @ 2:34
- recrystallized using 2-solvent method

(g) Dry weight c1 = 0.47g

c2 =

M.P. (°C) c1 = ~145

c2 =

% recovery c1 = $\frac{.47}{.91} = 51.6\%$

c2 = .91

(g) c1 + c2 =

60 mL added

Pre-Lab: Recrystallization Methods 22

Procedure

Data

Part III
cont.

Carry out M.P. Analysis

Approx M.P. $t_1 = 115^\circ$ $t_2 = 112.6^\circ$ DEM
 Actual M.P. $t_1 = 112.6^\circ$ DEM
 Literature M.P. = $113 - 114^\circ$

Carry out infrared spectroscopy


Notes

- 60ml H₂O added with no results of cloudiness...
- I dunked it in an ice bath to try and get crystals
- I attempted to bring out the crystals by boiling off MeOH

Determine unknown across ref from given M.P. & structure.

Unknown = - Reduced by half
 DEM ~ 25 mL - ~ 60 mL


- weight of final vial + cap = 17.98g
- weight of final v + cap + uk = 18.45
- Dunked reduced solution in ice bath after cloudiness formed
- Crystals formed @ 3:55
- performed vacuum filtration on the crystals
- vacuum filtration ended @ 4:20



Report: Recrystallization Methods

Table 1: Summary of Results

Identity	Acetanilide
M.P. °C	Observed: 112.6 Literature: 113
Color Before	A tan, sand colored granule
Color After	A colorless, fragmented crystal



$$\% \text{ yield} = \frac{.47}{.91} = 51.6\%$$

Possible Error during Experimentation/Conclusions

Part I: During the solubility testing, an error that could have been made was the misinterpretation of insoluble at R.T. and SS at R.T. since it was purely a qualitative observation. Many of the samples were small so this mistake could have been made.

Part II: A large source of error during this experiment was the fact that the filter in the Büchner funnel busted a hole during filtration and I lost yield. During this experiment I seemed to struggle a bit due to the amount of estimation and lack of precision during a recrystallization... since I had never done one before, it was challenging.

Part III: During this experiment, I added too much "poor" solvent and ended up boiling off a lot to reach a ~~cloud~~ cloud point in solution. I believe the amount of "boiling off" that was done affected my yield.

Conclusion of M.P. Analysis/spectroscopy ~~~~~>

Report: Recrystallization Methods


Part III conclusions cont.

M.P. Analysis: After obtaining the M.P. of 112.6°C , Acetanilide was the most likely contender as ~~my guess~~ the unknown, but benzoic acid could have also been possible if my purity was off due to increased impurities lowering the M.P. of a substance... to be 100% sure of what the unknown was, spectroscopy was performed.


Infrared Spectroscopy: Upon investigation of the graph, there were three main features.

- A spike at 1700cm^{-1} which indicates a carbonyl group
- High activity at $700-1000\text{cm}^{-1}$ which indicates a benzene ring
- A broad spike at 3400cm^{-1} which indicates an amine

After taking this information into consideration, the unknown was deduced as Acetanilide due to the amino group being present as well as both common groups shown above.



vs.



Acetanilide


Benzil

x David Miller, 10/19/22

Image 2: Comparison of Acetanilide and Benzil

[Return To Table of Contents](#)

EXTRACTION METHODS USING ACIDS AND BASES



Pre-Lab: Extraction Methods 23

• Objective

- To separate a solution of Benzil and 4-methylbenzoic acid using extraction methods and to recrystallize the respective compounds.

• Theory

- Extraction is an important technique that is used by chemists to separate out two different organic solutes in the same solution. This is done by using an acid or base to "attack" the protonation site of the molecule, thereby making it slightly polar \pm through dipole moments and therefore miscible in a polar solvent such as H₂O. This is important because it allows the chemist to separate the two solutes and perform various techniques on either one. (recrystallization/^{or} ~~reprecipitation~~ red protonation)
_{distillation}

- The most efficient way (for yield) is to perform multiple extractions... This can be demonstrated by the partition coefficient which, when calculated using multiple yields will read as follows...

$$K_p = \frac{\text{Amount A in Organic solvent} / V(\text{mL}) \text{ of organic solvent}}{\text{Amount A in Aqueous solvent} / V(\text{mL}) \text{ aqueous solvent}} = 20$$

Cont. \rightarrow

Pre-Lab: Extraction Methods

24

- If you enter in $A = 2.0g$ and compare it with a single & double extraction you get as follows...

Single


$$K_p = \frac{x/100(\text{ml})}{(2.0-x)/100(\text{ml})} = 20$$
$$x = 1.90$$

Double

$$K_p = \frac{x/50(\text{ml})}{(2.0-x)/100(\text{ml})} = 20$$

$$x_1 = 1.82g \quad x_2 = 0.16g$$

$$\underline{\underline{x_{\text{tot}} = 1.98g}}$$



- After this phase is complete and both compounds are back in their organic layers, a wash is performed to eliminate any traces of H_2O in the organic layer... Drying agents can also be used to eliminate unwanted aqueous layer.

Pre-Lab: Extraction Methods

25

Diagrams





Diagram 1: Separatory Funnel

Diagram 2: Gravity Filtration

Tables of Reagents & Products

Name	Structure	M.P.	Color	Hazards
Benzil	<chem>O=C(c1ccccc1)c2ccccc2</chem>	94-97°C	yellow (solid)	<ul style="list-style-type: none"> Causes serious eye irritation May cause resp. irrit.
p-Toluic acid	<chem>CC1=CC=C(C(=O)O)C=C1</chem>	215-217°C	White (solid)	<ul style="list-style-type: none"> serious eye/resp irrit. skin irrit. harmful if swallowed
HCl Hydrochloric acid	$H-Cl$ $H^+ \quad Cl^-$	-35°C	colorless (liquid)	<ul style="list-style-type: none"> Causes severe burns strong acid
Sodium bicarbonate	<chem>NaHCO3</chem> <chem>[Na+].[O-]C(=O)O</chem>	270°C	White (solid pow)	<ul style="list-style-type: none"> It's baking soda, you're good
Ethyl Acetate	<chem>CCOC(=O)C</chem>	-83.5°C	Colorless (liquid)	<ul style="list-style-type: none"> Flammable may cause dizziness eye irrit.
Water (l)	<chem>O</chem>	100°C	colorless liquid	<ul style="list-style-type: none"> may choke if you suck at drinking

Pre-Lab: Extraction Methods

26

References

~~Fishersci.com~~


1. Fisher Scientific (NID). Fishersci.com
MSDS (Table of reagents)
2. LabFlow.com (2022). Experiment 3, Extraction Methods
Storage-cdn, labflow.com

Definitions of Terms

Extraction: The moving of the desired compound from one layer to the other in immiscible solvents through the use of acids and bases.

Wash: The use of salts to remove impurities from the organic layer and transport them to the aqueous layer.

Procedure on next pg



Pre-Lab: ^{Dist} Extraction Methods

27

Table of Drying agents

Agent	Advantage	Disadvantage
Anhy MgSO ₄	Rapid high cap	Reacts w/ alcohols & amines
Anhy CaSO ₄	Very powerful	" " + low H ₂ O cap
Anhy Na ₂ SO ₄	Chemically inert	Slow, best for initial drying
Anhy K ₂ CO ₃	Rapid	Reacts w/ acids
Anhy CaCl ₂	Rapid	Reacts w/ alcohols

Procedure

Data

- Obtain 1.0g mixture of approx 0.5g each of 4-meth... & Benzoil weight
- Dissolve mixture in about 25ml of ethyl acetate
- 100ml beaker
- 125ml sep funnel
- powder funnel
- After dissolved, pour beaker contents into the 125ml sep funnel
- Rinse beaker w/ 5ml more acetate
- Rinse the powder funnel to add to sep funnel
- Extract using 10ml 125ml NaHCO₃ (saturated) and drain into a 250ml beaker
- repeat 2 times
- Chill the 250ml beaker

$$K_p = \frac{x / 30 \text{ ml}}{(.97 \text{ g} - x) / 10 \text{ ml}} = 20$$

$$\begin{aligned} \text{Vial} + \text{cap} + \text{powder (g)} &= 18.78 \text{ g} \\ \text{" " - powder} &= 17.81 \\ \text{powder} &= .97 \text{ g} \end{aligned}$$

- added 25ml Ethyl acetate + 5ml 30ml total
- Added NaHCO₃ to organic layer
- performed extraction and separated into a flask 2x
- Added HCl (10-15ml) to the soln. (aq)
- Added NaCl to the organic soln to wash the water out

Pre-Lab: Extraction Methods

28

Procedure

- Neutralize the chilled solution using 1M HCl. Stir vigorously.
- Do I see a precipitate?
- Is it slightly acidic to the pH paper? ICE IT
- While \uparrow is cooling, add about 10mL of saturated NaCl(aq) to the organic solution to remove H₂O.
 - Drain excess into a new beaker
 - Label beaker "NaCl soln"
- Transfer clean organic ~~phase~~ solution to a 125mL E. Flask and dry w/ MgSO₄.
 - Cork the flask
 - Swirl flask
 - Inspect MgSO₄
- Use vacuum filtration through a buchner funnel to obtain precipitate from chilled soln
 - Rinse w/ DI water (cold)
 - Make sure DI water is ICE
- Decant organic layer through gravity filtration (Diagram 2) and store in a 250mL RB flask.
- Discard drying agent!

Data

- Transferred Benzil soln. to the 125mL E-Flask
- Added drying agent after pulling H₂O out of the bottom of the flask
- After adding adequate HCl to aqueous soln. a lot of precipitate formed
- I removed precipitate using vacuum filtration and rinsed with icy DI water.
- I performed a gravity filtration on the organic solution.
- Cleaned up station
- Evaporated off the ethyl acetate in the roto vac device.

Pre-Lab: Extraction Methods

29

WEEK #2 Procedure

Data
Benzil

✓ Benzil

○ obtain a gross weight of RBF flask and its contents.

○ add 10 mL methanol and heat on a hot plate until residue dissolves.

○ Transfer solution to a 50 mL E-Flask.

○ Rinse w/ 5 mL MeOH to get all the benzil

○ Add a boiling stick and bring to a simmer on a hot plate.

○ perform a hot grav. filt. using a Kim wipe pipette (Monstener)

○ concentrate to about 10 mL and add MeOH until a cloud point is reached

○ add MeOH until soln. has cleared and allow to cool to RT.

○ If Benzil "oils out" of soln. add a bit more MeOH and repeat boiling process

○ Chill in an ice bath and vacuum filter crystals

○ Obtain MP

○ Obtain mass

Weight RBF = 135.70 g

RBF = 136.05 g

○ Gathered supplies

○ Started heating process 1:55 p

○ began hot grav. filtration 2:00

- Reduced filtrate until about 10 mL remain

- Added MeOH until cloud point was reached

- dispersed clouds w/ MeOH

- placed benzil solution in ice bath until crystals formed

- Used vacuum filtration to remove crystals and stored on watch glass.

- I lost a bit from the solvent used to "stick" the filter

4-Me Benzoic acid

- repeated process of recrystallization used previously

Notes that differ from previous

- The cloud point method didn't go as smoothly as benzil.

Pre-Lab: Extraction Methods

30

Procedure

Data

① 4-Methyl benzoic acid
Orecrystallize using
an EtO/MeOH two
solvent system.
O obtain MP. & mass

opresent both vials
to instructor.

- Performed MP analysis

4-Me = 180-181°C

Benzil = 94-97°C

M.P. Benzil = 95.8°C

Mass Benzil = .46g

% yield Benzil = 92%

M.P. 4-Me B-Acid = 187.1°C

Mass " " = 47.95 - 47.81 = 0.14g

% yield " " = 28%

Report: Extraction Methods

296

Summary Table



Name	4-Methylbenzoic acid	Benzil
Observed M.P.	181.7°C	95.8°C
Observed color	colorless	181.4°C ^{tan} yellow
% yield	28% ^{max}	92%
Structure		

Table 1: Summary Table for extraction Methods

Critiques

Benzil ~ The recrystallization of benzil went very well with a 92% yield... My only critique with this is the fact that the bichner funnel was not dampened and I lost yield into the waste solution.

4-Me benzoic acid ~ Upon extraction of the 4-Me benzoic acid yield must have been lost. I believe the most yield was lost when adding the HCl to obtain the precipitate... I either added too much or too little. I suspect something went wrong at this step because the recrystallization process went fairly smooth. Overall, these mistakes majorly hurt the yield making it 28%.


x. David Miller

11/1/22

[Return To Table of Contents](#)



OPTICAL RESOLUTION OF (+/-)MANDELIC ACID



Optical Resolution of (\pm) Mandelic Acid

31

Objective: To separate a racemic mixture of \pm Mandelic acid into its respective enantiomers using extraction, recrystallization, and rotary vacuum evaporation. The purpose of this is to obtain the optical resolution of the (R) and (S) enantiomers.

Theory: To meet the objective of separation of enantiomers, one must manipulate the chemicals in a way that lets them exploit a difference... This is done by using (+)-cinchonine to act as a base and deprotonate (R) and (S) Mandelic acid. By doing this, the enantiomers that once had the same properties become diastereomers with differing properties regarding solubility. This allows us to separate the diastereomers using recrystallization. After this, the diastereomers must be reprotonated to become (R) and (S) enantiomers again, so HCl is added to the mixture of the deprotonated mandelic acid still in solution with the cinchonine. After this, they are separated through extraction due to differing solubilities. After the (R) and (S) enriched solutions of Mandelic acid have been synthesized, an optical resolution is performed to calculate the specific rotation and % enantiomeric excess of the enantiomers. The reason optical resolution is so crucial is because of how enantiomers will bend light opposite of each other in either a dextrorotary or levorotary direction. Using this information, we can calculate how pure the enantiomer in a solution is.

Eqs, Diagrams, Mechanisms

Pre-Lab! Optical Resolution of Mandelic Acid

Theory cont.







IMAGE 1: Gravity Filtration



IMAGE 2: Separatory Funnel





Enantiomers: Same physical props.
 Deprotonation / Protonation



Diastereomers: Different phys properties... separate using recryst.



less soluble diastereomer



Pre-Lab!
Optical Resolution of (+) Mandelic Acid

33

Theory cont. : Equations

Theoretical Yield

Cinchonine = 1.5g
Mandelic Acid = (R) + (S) = .75g

Specific Rotation $[\alpha]_D = \frac{\alpha}{cl}$

α = Degrees rotated c = g substance / mL solvent l = 1 dm = 20 cm
(Mand Acid) (Et acetate)

% Enantiomeric Excess

$$\% EE = \frac{[\alpha]_D \text{ sample}}{[\alpha]_D \text{ from major Enant}} = \frac{\text{mol major E} - \text{mol minor E}}{\text{Total mol both}} \times 100\% = \frac{\% \text{ major} - \% \text{ minor}}{100\%}$$

Calculating $\% R$ and $\% S$

if... pure S = +50 and rotation mixture = -40 then $\% EE = \frac{-40}{-50}$
pure R = -50







\downarrow
 $\% = 80\%$ if $R + S = 100\%$ and $R - S = 80\%$ then $R = 100\% - S$

\downarrow
and therefore $100\% - S - S = 80\%$ $= -2S = -20\%$ $= \frac{-20\%}{-2}$ $S = 10\%$
 $R = 90\%$

TABLE OF REAGENTS NEXT \rightarrow

Optical Resolution of (+) Mandelic Acid

Table of Reagents

Name	Structure	MMg/mol	Hazards info
Water		18.02g	May cause a sudden urge to urinate
EtOH Absolute		46.08	<ul style="list-style-type: none"> Highly Flammable Causes serious eye irritation
Ethyl Acetate		88.11	<ul style="list-style-type: none"> Highly Flammable Causes serious eye irritation May cause drowsiness or dizziness
HCl	H-Cl	36.46	<ul style="list-style-type: none"> May be corrosive to metal Highly acidic Causes severe burns
Cinchonine		294.4	<ul style="list-style-type: none"> Harmful if swallowed May cause an allergic skin reaction.
Mandelic Acid (+)		152.15	<ul style="list-style-type: none"> Corrosive Causes serious eye damage
Magnesium Sulfate		120.36	Non Hazardous

References

- Labflow - Optical Resolution of Mandelic acid; Enantiomeric excess (2022). <https://www.Labflow.com/app/course/2063>
- Experiment 4 - Optical Resolution. Labflow.com (2022). <https://www.Labflow.com/app/course/2063>
- Fisher Scientific MSDS - H₂O; EtOH; HCl; MgSO₄; cinchonine; Mandelic acid; Ethyl Acetate. (2020). FisherSci.com

Procedure

Data

~~DEM~~

Pre-Lab
Optical Resolution of ¹³C Mandelic Acid

36

Procedure

- Obtain 0.75g racemic and 1.5g cinchonine from Dr. G and put into a 125 mL E flask.
- Add 30 mL Et Acetate
- Perform a 1 solvent recrystallization & completely dissolve at B.P
- Reduce until a small amt. of crystals form, ~40 mL left in ice bath!
- Vacuum filter crystals
- Repeat above process using ^{fresh} the dried crystals.
- Recrystallize in the 125 mL EA.
- Separate using extraction 2x
 - o organic = 2x 25 mL Et Acetate
 - o (aq) = 25 mL 3M HCl
 - o add HCl first! 25 mL
- Transfer organic layers to a 200 mL E flask
- Drain (aq) phase into a 250 mL beaker for both extractions
- Dry the organic phase w/ MgSO₄
- Perform gravity filtration
- Neutralize cinchonine w/ 3x 26 mL NaOH and vacuum filter the precipitate

Data

- Actual weight MA = 0.74g
- Actual cinchonine = 1.52g
- Dissolved in about 55 mL Et acetate
 - o 30 + 25 mL = 55 mL
- Reduced down to 30-40 mL (it was cloudy)... Recrystallized
- Performed vacuum filtration
- repeated recrystallization @ about 35-60 mL
- Crystals ~~formed~~ (salt) cinch + mand
 - enriched in (R)
- Measured out HCl (25 mL)
- added HCl to mand acid
- After extraction I reprecipitated the cinchonine with NaOH and set the (R) mand acid aside as the organic phase.
- I vacuum filtered cinchonine to get my yield after it was pH 7
- Weight RBF = 134.98g
- neutralized (aq) phase with NaOH... it was slightly basic but I proceeded anyways, ~28.5 mL
- Obtained precipitate w/ vacuum filter

Procedure

- Recrystallize the precipitate using a MeOH/H₂O scheme
- Vacuum filter crystals and weigh/turn it in with the report sheet.
- Use the Rotovac to evaporate the solvent to leave only the mandelic acid.
- Add 10ml Absolute EtOH to the flask to dissolve all contents.
- Transfer to a clean sample vial via a clean pipette.
 - if soln. is cloudy, filter through filt. paper.
- Obtain a polarimeter cell...
 - Rinse w/ a pipet of solution
 - Discard rinse
 - Fill to the brim like videos
- Obtain info for specific rotation
- Discard contents of cell
 - Organic waste
- Rinse cell w/ EtOH
- Calculate the following

Data

- performed recrystallization of cinchonine using MeOH/H₂O
 - Added 35ml Me to start
- Mass cinchonine = 0.4g
 - Added 15ml more
 - Added 30ml more
 - Reduced to about 40ml
- Mass Crude Mand Acid = 0.4g
- Concentration of Solution = 0.4g mandelic acid / 10ml EtOH
 - Added H₂O to recrystallize cinchonine
 - Performed vacuum filter (10ml)
 - Added EtOH to MA (R enriched)
 - Obtained observed rotation
- Disposed of MA and placed cinchonine in a vial
- Performed a M.P. Analysis on cinchonine (+)
 - MP = 248.8
 - it burned, but didn't melt haha.

Amount polarized = -1.45°
 Temp of room ($^\circ\text{C}$) = 21.4°C
 Pure (S) = $+160^\circ$ in EtOH Pure (R) = -160°
 Specific Rotation = -36.25°
 $\% \text{ EE} = 22.66\%$
 $\% \text{ (S)} = 36.7\%$
 $\% \text{ (R)} = 61.39\%$

Report: Optical Resolution of Mandelic Acid 36b

Critique of Experiment

- The experiment seemed to go smoothly for the most part. A perfect yield was not obtained so I will begin by coming up with possible explanations. There is a good chance I lost yield during recrystallization of mandelic acid since the solubilities of (R) and (S) are very similar when (R) is bound to cinchonine.
- My yield for cinchonine was around 40% which was not great. I believe there was error during both recrystallization and reprotonation to form the precipitate.

November 16th '23

Xavier Mills

Summary Table for (R) Mandelic Acid


Theoretical yield	.75g (R+S racemic) / 2 = .375g
Observed Rotation	-1.45°
Specific Rotation $[\alpha]_D^{25}$	-36.25°
Enantiomeric Excess	% EE = 22.66%
%(R) & %(S)	61.3% (R) 38.7% (S)



DEHYDRATION OF TETRAHYDROLI NALOOL.



(ALKANE \rightarrow ALKENE)




Pre-Lab: Dehydration

38

Objective: To use an acid to dehydrate a tertiary alkene alcohol to create and isolate a desired alkene.

Theory: In a dehydration, ~~one uses~~^{DEH}, an E1 mechanism occurs which is a multi-step process that "cleaves" the leaving group and forces the carbocation that is left to form a double bond through "kicking out" the β hydrogen. After these alkenes are formed, they are then separated using distillation since their boiling points have changed. Since tetrahydrofuran has 5 possible isomers (constitutional), we must determine the percentages of trisubstituted vs. Di".



Appendix 1: Distillation

Pre-Lab: Dehydration

Theory cont.



39

Theoretical yield

$$\text{MW tetrahydrothalool} = 158.285 \text{ g/mol}$$
$$\text{mol " " } = \text{g sample} \times \frac{1 \text{ mol}}{158.285 \text{ g}} = \text{mol sample}$$
$$= \text{Theoretical yield}$$


$$\frac{\text{mol product}}{\text{mol sample}} \times 100\% = \% \text{ yield.}$$

Dehydration of tetrahydrothalool ~ E1



Trisubstituted = Δ
Disubstituted = \sim
Entgegen = E
Zusammen = Z

5 Isomers



Mechanism 1: Isomers

Pre-Lab: Dehydration

40

Table of Reagents


Name	Structure	B.P. (°C)	Hazards
Sulfuric Acid	$\text{HO}-\overset{\text{O}}{\parallel}{\text{S}}-\text{OH}$	290	• Strong acid! • Cause's severe burns • Respiratory irritation.
Myrcene (similar to the products)		165°C	• Flammable • Skin irritant • may cause cancer • eye irrit
Magnesium ²⁺ Sulfate ²⁻	$\text{Mg}^{2+} \text{SO}_4^{2-}$	N/A	• Harmful if swallowed • may cause eye irritation

Table 1: Table of reagents

Procedure

Misc en place:

- 100 mL RBF
- 50 mL RBF
- Distillation setup
- Wood blocks
- Ice bath
- Boiling chips!
- 10 mL grad cyl
- Monstener pipette + Kim Filter
- pipettes
- vial
- Erlenmeyer (for drying)

Pre-Lab: Dehydration

41

Procedure

- Clamp a 100mL RBF in Hood station
- Take a 10mL GCyl and weigh out exactly 5.0mL in grams
- Transfer to RBF
- Clean all equipment used very well
- Measure out 10mL 6M H_2SO_4 ... use funnel to pour into the RBF
- Add boiling chips
- Set up the RBF just like "Apparatus #1" for distillation.
- Crank VARIAC to 90-100%, at 120V and dial back to a simmer after boiling occurs.
- RECORD temp at first drop distilled
- Distill until most of the liquid is gone.
- Remove aqueous layer in the 50mL collection RBF using a pipet and transfer to an E-flask...
 - Dry w/ $MgSO_4$

Data

- G-cyl mass = ~~4.24g~~
- Vial mass = ~~17.72g~~ ~~20.28g~~
- G-cyl + tetra = 49.11g
- tetra = ~~20.28g~~ 4.13g
- Product = ~~20.28g~~ 2.56g
- Set up RBF and distillation setup
- Measured out 5.0mL of tetrahyd
- DONT FORGET TO NEUT. THE ACID!
- Measured out 10mL H_2SO_4
- added to RBF
- Started distillation.
- Temp at first drop = 84.5°C
- Temp at stop of dist = N/A
- After removing H_2O , it is a bit soapy and forms bubbles.
- I dried the organic phase out
- filtered through a monster phette

Pre-Lab: Dehydration

42

Procedure

Data

○ obtain a tare weight of the sample vial and setup a constant gravity filtration.

○ obtain weight of product obtained and compare to T-yield to get % yield.

○ Rinse pipet with acetone before disposal.

○ Get an IR spectrum of the product

○ Get a reading from the NMR spectrometer

○ Turn in product.

- Obtained IR
- evaluated spectra
- There seems to be little to no -OH bonds in product... this is indicative of alkenes

- Obtained NMR?
- evaluated shielded and unshielded protons and data collected.

$$\%D = 42.9\%$$

$$\%T = 57.1\%$$


$$\%D = \left(\frac{I_{4.7}}{I_{4.7} + 2I_{5.1}} \right) \times 100$$

$$\%T = \left(\frac{I_{5.1}}{I_{5.1} + 5I_{4.7}} \right) \times 100$$

$$\% \text{ yield} = 62.00\% = \frac{2.56 \text{ g product}}{4.13 \text{ g reagent}}$$

% Transmittance

100 DM Alkene Spectra



lack of
OH

← C-H
stretch


← C=C
Carbenes

← C=C saturation

Wavenumbers (cm⁻¹)

4000 3500 3000 2500 2000 1500 1000

2022-11-17



$$\%D = \frac{1.5}{(1.5 + 2(1.0))} = .429 = 42.9\%$$

$$\%T = 100\% - \%D = 57.1\%$$

-OH being low in IR
 Parallel the evidence
 of lower levels of
 alkanes.

David Miller

Acquisition Parameters

- Acquisition Time: 2022/11/17 13:33:18
 - Nucleus: ¹H
 - Solvent: Chloroform-d
 - TX Frequency: 60.53 MHz
 - Scans: 16
 - Time Per Scan: 4.70 s
 - Spectral Width: 12 ppm
 - Dwell Time: 0.29 us
 - Digital Resolution: 0.04 Hz
 - NPoints (Complex): 2048
 - Zero Filling: 7.00
 - Apodization: 0.10
 - Receiver Gain: 40.89
 - Pulse Width: 14.4 us
- Peaks:
- 1: -63.5 Hz (-1.05 ppm)
 - 2: -62.4 Hz (-1.03 ppm)
 - 3: -60.9 Hz (-1.01 ppm)
 - 4: -59.5 Hz (-0.98 ppm)
 - 5: -57.2 Hz (-0.95 ppm)
 - 6: -55.7 Hz (-0.92 ppm)
 - 7: -51.6 Hz (-0.85 ppm)
 - 8: -48.7 Hz (-0.80 ppm)
 - 9: -45.6 Hz (-0.75 ppm)
 - 10: -43.1 Hz (-0.71 ppm)
 - 11: -42.0 Hz (-0.69 ppm)
 - 12: -35.4 Hz (-0.59 ppm)
 - 13: -1.0 Hz (-0.02 ppm)
 - 14: 1.4 Hz (0.02 ppm)
 - 15: 3.3 Hz (0.05 ppm)
 - 16: 49.0 Hz (0.81 ppm)
 - 17: 54.8 Hz (0.90 ppm)
 - 18: 66.3 Hz (1.10 ppm)
 - 19: 68.2 Hz (1.13 ppm)
 - 20: 72.6 Hz (1.20 ppm)
 - 21: 75.2 Hz (1.24 ppm)
 - 22: 78.9 Hz (1.30 ppm)
 - 23: 80.8 Hz (1.34 ppm)
 - 24: 91.0 Hz (1.50 ppm)
 - 25: 95.3 Hz (1.57 ppm)
 - 26: 104.5 Hz (1.73 ppm)
 - 27: 107.5 Hz (1.78 ppm)
 - 28: 109.0 Hz (1.80 ppm)

1A Report: Dehydration of Alkanes 41b

Summary

Procedure: My yield was not perfect, but it was decent. The only thing that I believe could have hurt my yield was the drying process as I added quite a bit of $MgSO_4$.

IR spectroscopy: Upon observation of the IR, the experiment seems to be successful. There was a lack of an "OH bulge" at $3200-3500\text{ cm}^{-1}$ which indicates dehydration was successful. There was also significant evidence of $C=C$ due to a reading at around 1650 cm^{-1} .

NMR: The NMR confirms the success due to the deshielding downfield on the graph which indicates the presence of alkenes with low levels of alkanes.

Summary Table

A-yield	2.56g
T-yield	4.13g
% yield	62%
%D	42.9%
%T	57.1%
Observed BP	84.5°C

X David Miller

11/30/2022


[Return To Table of Contents](#)



WILLIAMSON ETHER SYNTHESIS



THE ETHER BUNNY



Pre-Lab: Will Ether Synthesis

43

Objective: To utilize a stoichiometric mixture of eugenol, NaOH, and ethylene glycol ditosylate to create an ether product through refluxing the reaction.

Theory:



The following reaction takes place in a reflux of reactions which is when a constant flow of energy is provided to the reaction via boiling, but none of the vapors are boiled off; instead, they are condensed back down into the flask. Today's experiment will be carried out in a safer manner than NaCS being added by using NaOH to form a salt out of the eugenol. This salt is then used to "cleave" the two halves of the ditosylate to create the Ethylene glycol Ether. This happens due to the two halves being held together are "tosyl" groups which are good leaving groups for an S_N2 reaction to occur.

To analyze the samples, a nuclear magnetic resonance will be taken which is a spectroscopic procedure in which the frequencies of the molecules resonating is measured and graphically presented.

Pre-Lab! Williamson Synthesis of Ethers


44

Equations/Mechanisms




Equation 1: Basic Williamson Ether synthesis using Na metal.

A way to get around using Na(s) is to do the following reaction (much safer)



OT = Tosyl group



Good leaving group

Pre-Lab: Williamson Ether Synthesis

45

Equations needed (Theory)

$$\% \text{ yield} = \frac{\text{actual yield}}{\text{Theoretical yield}} \times 100\%$$

Calculations for stoichiometry

$$g \text{ diethylate} \times \frac{\text{mol diethylate}}{g \text{ diethylate}} = X \text{ mol}$$

$$3X \text{ mol} \times \frac{g \text{ Enguol}}{\text{mol Enguol}} \times \frac{\text{mL Eng}}{g \text{ Eng}} = \text{mL needed}$$

$$3 \cdot 3X \text{ mol NaOH} \times \frac{g \text{ NaOH}}{\text{mol NaOH}} \times \frac{\text{mL NaOH}}{g \text{ NaOH}} = \text{mL needed}$$

Freq. independent & scale Formula

$$\delta \text{ ppm} = \frac{\nu \text{ sample (Hz)} - \nu \text{ reference (Hz)}}{\nu \text{ spectrometer (MHz)}} \times 10^6$$

Terms:

Upfield: Right part of spectrum... contains "shielded" protons.

Downfield: Left part of the spectrum... contains "deshielded" protons.

Diagram for reflux




Figure I: Reflux setup

Pre-Lab: William Ether Synthesis

46

Table of Reagents/Products





Name	Structure	M.P./BP	M.W.	Hazards
Eugenol		-11°C / 254°C	164.2 g/mol	<ul style="list-style-type: none"> Harmful if swallowed May cause allergic rxn Causes eye irrit
Sodium Hydroxide	$H-O^- Na^+$	450°C / 2000°C	40 g/mol	<ul style="list-style-type: none"> Corrosive to metals Causes severe burns Resp. Irritation
Ethanol		20°C / 77.1°C	46 g/mol	<ul style="list-style-type: none"> Flammable Suspected of causing cancer Damages fertility
Ethylene glycol bis-p-toluenesulfonate		122°C / N/A	310.45 g/mol	<ul style="list-style-type: none"> Skin irritation eye " respiratory "
Methanol	$OH-CH_3$	-98°C / 64.7°C	32.04 g/mol	<ul style="list-style-type: none"> Dang to organs Resp. irrit Flammable Toxic if swallowed
Ethyl Acetate		-85°C / 76°C	88.11 g/mol	<ul style="list-style-type: none"> Flammable May cause dizziness eye irrit
X	X	X	X	X

Table 1: Table of reagents

Pre-Lab: Williamson Ether Synthesis
Procedure

47

Data

Measure out 1.25 mL (1.378g) of eugenol into a 10 mL G-cyl
use a pipet

put eugenol in RBF w/
stir bar

Add .356 g (.008904 mol) NaOH to the RBF

Rinse into the flask with
4 mL EtOH

turn on magnetic stirrer
stir for 10 min @ RT

observe yellow color

obtain approx 1g dithosylate
Add to RBF

Add 4 mL EtOH to rxn

Reflux for 30 min

note exact time reflux
begins (boiling start)

pour mixture while still hot
into a sep funnel w/ 25 mL H₂O
remove stir bar

Rinse RBF using H₂O & ethyl
acetate into the sep funnel

extract w/ 25 mL Et acetate

Measured out .29g NaOH
Measured about 1.3 mL Eugenol

I began stirring at 1:32p
and noticed the orangy eugenol
change to a dark green.
- eugenol smells like cloves! yum!

added EtOH and dithosylate
started VARIAC on 40

Mass vial+cap+dtho = 18.58

Mass vial+cap = 17.60

Mass dithosylate = .98g

Reflux start: 2:11p

Reflux end: 2:45p

- Reflux mixture got sludgy...
I added 5 mL more EtOH.

- ended reflux

Next

Pre-Lab: Williamson Ether Synthesis
Procedure

48

Date

- Drain the (aq) phase
- Rinse organic phase with 15 mL Saturated (aq) NaCl
- Discard both (aq) phases
- Dry organic phase using $MgSO_4$
- Gravity filter out $MgSO_4$
- Evaporate filtrate using rotovap
- Recrystallize from methanol using a rocit/hot grav filtration (if needed)
- Perform vacuum filtration
- Obtain NMR spectrum
- Obtain M.P. of crystals.

- While waiting for the reflux to finish, I measured out the H_2O , NaCl, and EtOAc...
- I also measured out $MgSO_4$
- Folded a fluted filter paper
- Reflux started to get gooky and a very dark green.
- Reflux has a film on top
- Performed H_2O , EtOH, EtOAc extraction
- Dried organic phase
- disposed of (aq) phase
- Filtered organic phase
- rotovapped organic phase
- performed recrystallization in MeOH
- added .04g rocit
- treated again w/ rocit using .92g this time.
- After the hot grav filtration, I recrystallized the product.
- I vacuum filtered product for 15 min
- Crystals oxidized upon start of next lab
- Mass Crystals: 0.33g
- Di Ether M.P. = 82.5°C

FA Report: Will Ether Synthesis

47b

Critique: During the experiment, a few things occurred that I believe significantly affected my yield.

1. When doing the nitrite treatment, I believe too much was added and it removed some of my product.

2. During reflux, a reflux ring was not able to be observed which possibly means some of the product was boiled off and lost.

3. One last point is that the MP measurement might not be as accurate as it could be due to the thickness and static electricity of the crystals preventing them from properly falling into the MP analysis rod.

Summary Table

Mass Di Ether	0.33g
% yield	23%
MP. Di Ether	82.5°C
Color	Pink and white crystals.

X. David Miss

END OF ORGANIC CHEMISTRY LAB 2022

